

**3M™ CONCISE™ RESTORATIVE MATERIAL**

**Technical Profile**

# 3M™ CONCISE™ TECHNICAL PROFILE

<b><u>PRODUCT DESCRIPTION</u></b>	<b>3</b>
<b><u>COMPOSITION AND CURING MECHANISM</u></b>	<b>3</b>
<b><u>EASY HANDLING</u></b>	<b>4</b>
FINISHING AND POLISHING	4
<b><u>ASSETS</u></b>	<b>4</b>
COLOR STABLE	4
ADEQUATE SMOOTH SURFACE	4
MINIMAL STRESS IN RESTORATION	5
DENSITY	5
WATER SORPTION AND SOLUBILITY	6
COMPRESSIVE STRENGTH	6
<b><u>CLINICAL PERFORMANCE</u></b>	<b>6</b>
<b><u>MAINTENANCE AND REPAIR</u></b>	<b>6</b>
<b><u>REFERENCES</u></b>	<b>7</b>
<b><u>APPENDIX A</u></b>	<b>9</b>
<b><u>APPENDIX C. FINISHING AND POLISHING TOOLS</u></b>	<b>10</b>
<b><u>APPENDIX C. SUMMARY COMPOSITION AND PHYSICAL PROPERTIES</u></b>	<b>11</b>

# 3M™ CONCISE™ TECHNICAL PROFILE

## Product Description

Concise™ Restorative composite consists of two pastes which are mixed in equal volume proportions. Both pastes have the same concentrations of matrix resins and filler powder. The pastes differ from each other in that Paste A contains the accelerator and Paste B the catalyst.

Concise™ is a chemically activated, strong, esthetic and long-lasting anterior restorative material. Retention and sealing are provided through the use of the acid etch technique and an adhesive system.

Concise™ is recommended for use in Class III and V restorations, for limited use in Class I restorations in premolars not subjected to occlusal stress, and Class IV restorations. Concise™ contains 78% ( by weight) quartz filler with an average particle size of 9 µm ( the particle sizes ranges from 1 to 40 µm).

## Composition and Curing Mechanism

The major component is a dimethacrylate (Bis-GMA) described in U.S.Patent 3,066,112 by R.L.Bowen. The diluent resin used is also a dimethacrylate, triethyleneglycol dimethacrylate (TEGDMA). No methyl methacrylate is present. The resin part is about 22 percent by weight of the composite.

The filler powder is a form of ground quartz with average particle diameter of about 9 µm. A small percentage of microfill particles is added to the filler fraction (Lutz & Phillips, 1983; Pilliar, Vowles & Williams, 1987). It is treated with a silane so that it is chemically bound to the resin during cure. The solid fraction is about 78 percent by

weight of the composite. This corresponds with approximately 67 percent by volume (Willems *et al*, 1993; Cross, Douglas & Fields, 1983).

Other ingredients are present in very minor amounts. The catalyst is benzoyl peroxide. The accelerator is a tertiary amine of the dimethyl- paratoluidine type. Inhibitors are present in parts-per-million. Titanium dioxide and stable iron oxides are present for pigmentation and provide a consistent shade.

Approximately equal volumes of pastes A and B of Concise™ are mixed thoroughly for 20 seconds. The resultant paste has a working time of 60-90 seconds ( depending on temperature) after which time polymerization will have proceeded sufficiently far to turn the fluid paste to a gel. In two minutes the paste will have hardened sufficiently to allow removal of excess material, and after a further two minutes final finishing and polishing procedures can be carried out.

Simply described, when the two pastes are mixed, the accelerator immediately begins "activating" the catalyst causing it to form free radicals. The reaction of these radicals with the dimethacrylate resin molecules is temporarily prevented ( working time) because the radicals react first with the inhibitors. When the small amounts of inhibitors are consumed, the catalyst free radicals then react with the resins and the vinyl coating on the filler particles to form chains of resin and coated filler with cross-linking between chains beginning (gel time). This reaction continues during the hardening time, linking filler particles and resin molecules into a highly cross-linked system ready for finishing.

## Easy Handling

### ***Finishing and polishing***

Delay contouring and finishing for a minimum of five minutes after matrix removal. Contour and polish restorations using Sof-Lex™ discs and strip systems, tungsten carbide or diamond finishing burs and mounted points.

Albers (1992) extensively describes the different finishing and polishing tools and techniques. His tips give useful guidance in addition to the instructions for use of the Sof-Lex system. Albers' tips are summarized in the next paragraphs.

Slightly overfill the composite. This allows for removal of the outer air inhibited layer which is softer than the inner parts of the restoration. Use micron diamonds (40-60 µm) for bulk reduction on surfaces unreachable with discs. Use slow or stall out speeds with copious amounts of water. This provides a smooth surface with minimal resin damage. Diamonds are generally better than discs for placing surface texture.

Medium or course diamonds leave a rough surface. This can extend the following finishing and polishing time. Avoid stones when approaching the final contour because loosened particles (caused by heat and vibration) may result in more surface porosity.

Burs cut the composite surface, which increases the likelihood of resin matrix fatigue fracture.

Sof-Lex discs and strips give an excellent finish and polish. Use the coarse discs with water and a very light touch because heat and friction can weaken the composite and enamel-resin interface. Keep discs constantly moving on

the surface to prevent heat and flat spots.

Flexible discs cut composite more rapidly than enamel. Therefore, they can ditch composite very easily. Use discs with more rigid backing (Sof-Lex XT) when polishing margins of materials with different cutting rates (enamel/composite; metal/composite). Polishing pastes, flour of pumice, tin oxide, rubber wheels and 12-fluted burs generally increase the surface roughness.

Farah & Powers (1998) give a good update of currently available diamonds, carbides, discs, mounted points and pastes. Appendix B gives an overview of the available instruments and their respective areas of use.

## Assets

### ***Color Stable***

The color stability of Concise™ restorative to UV light satisfies the relevant ISO standard.

Contrary to light-cured composites, Concise™ restorative does not show perceptible color changes upon curing (Seghi, Gritz & Klim, 1990).

### ***Adequate Smooth Surface***

Generally, the scientific community relates a minimized plaque retention to a smooth surface. Benderli *et al* (1997) found that more plaque retention was observed on a proprietary microfill composite resin that releases fluoride compared to Concise™!

A smooth surface should remain smooth after functioning under oral conditions for some time. The 'golden standard' is the surface roughness enamel after occlusal enamel-to-enamel contacts. Willems *et al* (1993) report

# 3M™ CONCISE™ TECHNICAL PROFILE

for the intrinsic roughness of enamel a value of 0.64  $\mu\text{m}$ . Concise™ has an intrinsic surface roughness after tooth brush abrasion of 1.44  $\mu\text{m}$  (Willems *et al*, 1993), which indicates that it makes sense to re-polish the restorations at recalls (see section Maintenance and Repair).

A slightly rougher surface is a concern to some authors in that a restoration may unacceptably abrade an opposing tooth. Enamel against enamel represents the golden' standard. Concise results in the same tooth profile reduction as enamel (Ratledge, Smith & Wilson, 1994). Suffice it to say that a microfill product such as Silux Plus actually is better than enamel.

## **Minimal Stress in Restoration**

Stress built-up during polymerization of a composite restoration depends on several factors: An effective bond of the adhesive (a loose restorations has no internal stress), polymerization shrinkage, the modulus of elasticity as well as the water absorption and the related hygroscopic expansion of the composite. The polymerization shrinkage of proprietary contemporary composite materials is essentially in the same ballpark. The volumetric polymerization shrinkage for Concise™ is in the range of 1.8-2.4% (Matschinske, Tappe & Sandner, 1989; Matschinske, Tape & Matschinske, 1991), which is somewhat smaller than that of major part of the market products. At last, the mode of curing is important. Chemically curing materials such as Concise™ have lower internal stress levels developing than light curing materials of similar composition. Due to their lower curing rate, both the adhesive bond strength has more time to develop to adequate levels and internal stresses developing in the restorations can relax for a longer time.

## **Density**

The density (Figure 1) of a material represents the number of grams per  $\text{cm}^3$  (more generally expressed: mass per unit of volume).




Figure 1. Density (in g/cc) of Concise™ (C) and minimum (Min) and maximum (Max) of competitive traditional anterior and universal composite resin restoratives.<sup>1</sup>

When the price per gram is multiplied with the density, one obtains the price per  $\text{cm}^3$ . A cavity represents a certain volume that needs to be restored. So, to compare prices, it makes sense to compare them on a price/volume-basis.

Figure 2 gives a conversion factor that relates to the price per volume. When one defines the price of Concise™ as 1 (per unit of volume), one obtains the price per unit of volume by multiplying the price per gram with the conversion factor. For instance, in case of a product with a conversion factor 1.2, one needs to multiply the gram-price by 1.2 to obtain a comparable price on a cc(ml)-basis. Suffice it to say that Silux Plus is much lower priced than the price per gram suggests.

<sup>1</sup> Competitive data competitive anterior and universal products from Momoi & McCabe (1994) and De Gee, Feilzer & Davidson (1993)

# 3M™ CONCISE™ TECHNICAL PROFILE




Figure 2. Conversion from gram- to ml-price for of Concise™ (Con) and minimum (Min) and maximum (Max) of competitive composite resin restoratives.

## Water Sorption and Solubility

Water sorption and solubility are within the limits required by the relevant ISO-standard. Usually, market products are within the ISO acceptance limits. A negative clinical influence of a somewhat higher water sorption or solubility within the ISO acceptance limits is not known.

## Compressive Strength




Figure 3. Compressive strength of Concise™ (Con) and minimum (Min) and maximum (Max) of competitive traditional composite resin restoratives. Data from Willems et al (1993).<sup>2</sup>

Enamel and dentine have compressive strength numbers of 384 and 297 MPa respectively (Craig, 1989). Phillips (1982) reports 305 and 248 MPa as

<sup>2</sup> Internal data for Concise by 3M are higher (see Appendix C). This is probably due to test differences.

fracture strengths of a molar and premolar respectively. Clinical performance results confirm that the compressive strength of Concise™ is sufficient for its indication for anterior use.

## Clinical Performance

Concise™ is one of the best documented composite systems on the market. Literature gives numerous results of clinical results. The table in Appendix A gives a brief summary of published clinical reports.

The clinical performance not only depends on Concise™. Certain aspects strongly depend on the choice of the dental adhesive or selection of an elastic resin modified glass ionomer liner. For more information on these products, we refer to dedicated 3M product profiles.

## Maintenance and Repair

Esthetic restorations need to be maintained. Both by the dentist and the patient. Plaque should be routinely removed so as to avoid adjacent soft tissue becoming inflamed. Discoloration, staining, soft tissue response, fracture or other deterioration effects may impair the functional aspect of the restoration. Both dentist (craftsmanship, maintenance during check-ups) and patient (nutrition, oral hygiene, life style) play a pertinent role and may cause the restoration to fail after a short period of time. For instance, an undercured composite in combination with a specific diet (red wine, coffee, cola, tea) may give unacceptable discoloration. Alcohol consumption may soften the resin matrix of (undercured) composite or compromise a bonding agent. High alcohol mouth-rinses will have a similar effect and may adversely affect the esthetics. Periodic refinishing or polishing of the restoration surface by the dentist

# 3M™ CONCISE™ TECHNICAL PROFILE

will exert a positive influence. Proper brushing, flossing, oral irrigation and interdental stimulation form other useful tools of the 'maintenance armamentarium'.

Composite repair is a viable solution to correct surface discoloration of existing restorations, to restore small areas of (recurrent) caries along the margin of an otherwise sound composite restoration or in a situation when complete removal of a very large composite restoration would jeopardize the integrity of a tooth. Pounder, Gregory & Powers (1987) describe a viable repair technique for Concise™. They apply an intermediate enamel or dental adhesive bonding agent to a 600-grit rough surface.

The following paragraphs give some practical tips by Albers (1992).

Avoid mechanical cleaning devices on composites because they may badly "pit" the restoration. When there is a need for re-polishing, avoid coarse prophypastes.

Fillers of macrofilled composite resin restorative may be dissolved and pitted by acidulated phosphate fluoride (APF) gels. It would be best to use non-APF fluorides on patients with Concise™ restorations.

Where small repairs are feasible, one should be careful with large repairs. Christensen (1995) explains the problems very well: "After a few hours in the mouth, resin restorations cannot be expected to bond well to new repair resin, because the chemical activity of the original restorative resin has been exhausted. If resin repair is being considered, and either the repair segment or the remnant of remaining resin is not expected to be self-retentive within the

tooth, the entire restoration should be replaced. If the old restoration can be considered retentive by itself, irrespective of the new portion, a repair can be considered. Usually, total replacement is better than repair." For instance, Mitsaki-Matsou *et al* (1991) report indeed that the repair bond strength ranges from 19% to 52% of the strengths of the unrepaired resins after storing the samples at 37 °C for periods up to one year.

## References

- Albers HF. Finishing direct restorations. *Adept* 1992; **3**: 1-16.
- Baratieri LN, Monteiro S, Andrado MAC de. The "sandwich" technique as a base for reattachment of dental fragments. *Quint Int* 1991; **22**: 81-85.
- Benderli Y, Ulukapi H, Balkanli O, Külekçi G. In vitro plaque formation on some dental filling materials. *J Oral Rehab* 1997; **24**: 80-83.
- Christensen GJ. Resin restorations for anterior teeth. *J Amer Dent Assoc* 1995; **126**: 1427-1428.
- Cross M, Douglas WH, Fields RP. The relationship between filler loading and particle size distribution in composite resin technology. *J Dent Res* 1983; **62**: 850-852.
- Crumpler DC, Heymann HO, Shugars DA, Bayne SC, Leinfelder KF. Five year clinical investigation of one conventional composite and three microfilled resins in anterior teeth. *Dent Mater* 1988; **4**: 217-222.
- Davis RD. A clinical comparison of three anterior restorative resins at 3 years. *J Amer Dent Assoc* 1986; **112**: 659-663.
- Dijken JWV van. The effect of cavity pre-treatment procedures on dentin bonding: A 4 year clinical evaluation. *J Prosthet Dent* 1990; **64**: 148-152.
- Eames WB, Rogers LB. Porcelain repairs: retention after one year. *Operat Dent* 1979; **4**: 75-77
- Farah J, Powers JM. Composite finishing and polishing. *The Dental Advisor* 1998; **15**: 1-5.
- Faunce FR, Myers DR. Laminate veneer restoration of permanent incisors. *J Am Dent Assoc* 1976; **93**: 790-792.
- Hansen EK. In vivo cusp fracture of endodontically treated premolars restored with MOD

# 3M™ CONCISE™ TECHNICAL PROFILE

- amalgam or MOD resin fillings. *Dent Mater* 1988; **4**: 169-173.
- Hurwitz H. The mutilated dentition: an alternative approach. *J Amer Dent Assoc* 1980; **100**: 703-705.
- Jokstad A, Mjör IA, Nilner K, Kaping S. Clinical performance of 3 anterior restorative materials over 10 years. *Quint Int* 1994; **25**: 101-108.
- Käyser AF, Mentink AGB, Meeuwissen R, Leempoel PJB. Het klinisch gedrag van metaal-composiet stiftopbouwen. *Ned Tijdschr Tandheelk* 1992; **99**: 401-403.
- Lutz F, Phillips RW. A classification and evaluation of composite resin systems. *J Prosthet Dent* 1983; **50**: 480-488.
- Matschinke U, Tappe A, Sandner B. Zur Polymerisationsschrumpfung von Kompositen. *Zahn- Mund- Kieferheilk* 1989; **77**: 783-789.
- Matschinske U, Tappe A, Matschinske F. Die Messung der Polymerisationsschrumpfung von Kompositen mit einer neuen Methode. *Z Stomatol* 1991; **88**: 15-21.
- Mettler P, Friedrich U, Roulet J-F. Studie über die Abrasion von Amalgam und Kompositen im Seitenzahnbereich. *SSO Schweiz Monatschr Zahnheilk* 1978; **88**: 324-344.
- Mitsaki-Matsou H, Karanika-Kouma A, Papadoyiannis Y, Theodoridou-Pahine S. An in vitro study of the tensile strength of composite resins repaired with the same or another composite resin. *Quint Int* 1991; **22**: 475-481.
- Phillips RW. *Skinner's Science of Dental Materials*, ed 8. Philadelphia: Saunders, 1982.
- Pilliar RM, Vowles R, Williams DF. The effect of environmental aging on the fracture toughness of dental composites. *J Dent Res* 1987; **66**: 722-726.
- Pounder B, Gregory WA, Powers JM. Bond strengths of repaired composite resins. *Operat Dent* 1987; **12**: 127-131.
- Raadal M. Follow-up study of sealing and filling with composite resins in the prevention of occlusal caries. *Comm Dent Oral Epidemiol* 1978; **6**: 176-180.
- Roulet JF, Mettler P, Friedrich U Die Abrasion von Amalgam und Kompositen im Seitenzahnbereich - Resultate nach 3 Jahren. *Dtsch Zahnärztl Z* 1978; **33**: 206-209.
- Roulet JF, Mettler P, Friedrich U Studie über die Abrasion von Kompositen im Seitenzahnbereich - Resultate nach 3 Jahren. *Dtsch Zahnärztl Z* 1980; **35**: 493-497.
- Seghi RR, Gritz MD, Kim J. Colorimetric changes in composites, resulting from visible-light-initiated polymerization. *Dent Mater* 1990; **6**: 133-137.
- Simonsen RJ, Stallard RE. Sealant-restorations utilizing a luted filled composite resin: One year results. *Quint Int* 1977; **8**: 77-84.
- Simonsen RJ. Preventive resin restorations: three-year results. *J Amer Dent Assoc* 1980; **100**: 535-539.
- Simonsen RJ. In: Waning A. Do those restorations look as nice after 5, 10 or 15 years? *Fenestra* 1997; **8**: 58.
- Smales RJ, Creaven PJ. Evaluation of clinical methods for assessing the surface roughness of restorations. *J Prosthet Dent* 1979; **42**: 45-52.
- Smales RJ. Incisal angle adhesive resins: A 5-year clinical survey of two materials. *J Oral Rehab* 1983; **10**: 19-24.
- Smales RJ. Long-term deterioration of composite resin and amalgam restorations. *Operat Dent* 1991; **16**: 202-209.
- Smales RJ, Gerke DC. Clinical evaluation of four anterior composite resins over five years. *Dent Mater* 1992; **8**: 246-251.
- Smales RJ, Webster DA, Leppard PI. Predictions of restorations deterioration. *J Dent* 1992; **20**: 215-220.
- Smales RJ. Effects of enamel-bonding, Type of restoration, patient age and operator on the longevity of an anterior composite resin. *Amer J Dent* 1991; **4**: 130-133.
- Varpio M. Proximoclusal composite restorations in primary molars: A six-year follow-up. *ASDC J Dent Child* 1985; **52**: 435-440.
- Willems G, Lambrechts, Braem M, Vanherle G. Composite resins in the 21st century. *Quint Int* 1993; **24**: 641-658.

# 3M™ CONCISE™ TECHNICAL PROFILE

## Appendix A

Reference	Remark
Baratieri, Moneiro & Andrada (1991)	Case, reattachment fragments, 3 years
Crumpler <i>et al</i> (1988)	Class III, five years
Davis (1986)	Anterior, three years
Dijken (1990)	Abrasion/erosion, four years
Eames & Rogers (1979)	Porcelain repairs, one year
Faunce & Myers (1976)	Laminate veneer, case/technique
Hansen (1988)	Class II, fifteen years
Hurwitz (1980)	Interim post/composite, case/technique
Jokstad <i>et al</i> (1994)	Class III, IV and V, ten years
Käyser <i>et al</i> (1992)	Metal/composite build-up
Mettler, Friedrich & Roulet (1978)	Posterior, wear, one year
Raadal (1978)	Preventive restoration, 2.5 years
Roulet, Mettler & Friedrich (1978)	Posterior, wear, one year
Roulet, Mettler & Friedrich (1980)	Posterior, wear, three years
Simonsen & Stallard (1977)	Preventive restoration, one year
Simonsen (1980)	Preventive restoration, three years
Simonsen (1997)	Closed diastema case, nineteen years
Smales & Gerke (1992)	Class III, IV and V five years
Smales (1983)	Incisal angle restorations, five years
Smales (1991)	Anterior, sixteen years
Smales (1991)	Anterior, sixteen years
Smales, Webster & Leppard (1992)	Anterior, eighteen years
Varpio (1985)	Class II, pedodontics, six years

# 3M™ CONCISE™ TECHNICAL PROFILE

## Appendix C. Finishing and Polishing Tools

Available finishing and polishing armamentarium and its area of use (adapted from Farah & Powers, 1998)

<b>Tool</b>	<b>Shape or Abrasive</b>	<b>Area of Use</b>
<i>Gross Reduction</i>		
25-45 µm diamonds	Flame	Labial/Lingual
	Needle, needle	Labial/Proximal
	Needle with rounded end	Occlusal
8-12 µm fluted carbides	Concave, tapered	Cervical, labial
	Egg	Lingual
	Needle	Occlusal/Proximal
Discs	Aluminum Oxide	Labial/Lingual/Proximal
Strips	Aluminum Oxide	Proximal
	Diamond - 45 µm	Proximal
<i>Finishing/Polishing</i>		
8-25 µm diamonds	Flame, flame	Labial/Lingual
	Needle with rounded end, needle	Occlusal
	Needle	Proximal
30-fluted carbides	Egg	Lingual
	Needle	Proximal/Occlusal
Discs	Aluminum Oxide	Labial/Lingual/Proximal
Strips	Aluminum Oxide	Proximal
	Diamond - 15 µm	Proximal
Mounted points	Aluminum Oxide	Labial/Lingual/Occlusal
Paste	Aluminum Oxide	Where accessible

# 3M™ CONCISE™ TECHNICAL PROFILE

## Appendix C. Summary composition and physical properties

Aspect	Unit	Time	Concise™
<b>Resin</b>			Bis-GMA TEGDMA
<b>Filler</b>			
<i>Type</i>			Silica (Quartz)
<i>Loading</i>	Vol. %		62
	Wt. %		78
<b>Density</b>	g/cc		1.86
<b>Compressive Strength</b>	MPa	1 h	260
<b>Diametral Tensile Strength</b>	MPa	24 h	55
<b>Water Absorption</b>	μg/mm <sup>3</sup>		<50
<b>Visual Opacity</b>			0.35
<b>UV Color Stability</b>		24 h	Pass
<b>Radiopacity</b>			No