
3M ESPE

Elipar™ FreeLight™ 2

LED Curing Light

Lampe de photopolymérisation LED

Luz LED para polimerización

Operating Instructions
Notice d'utilisation
Instrucciones de manejo

04/07
3M ESPE AG
Dental Products

Elipar™ FreeLight™ 2

LED Curing Light

Safety

PLEASE NOTE

Prior to installation and start-up of the unit, carefully read the instructions provided herein.

As with all technical devices, the proper function and safe operation of this unit depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.

1. The unit must be used in strict accordance with the following Operating Instructions. The manufacturer accepts no liability for any damage resulting from the use of this unit for any other purpose.
2. Prior to start-up of the unit make sure that the operating voltage of the Elipar FreeLight 2 (stated on the rating plate) is compatible with standard voltage. Operation of the unit at a different voltage may damage the unit.
3. Do not look directly at the light guide and do not stare at light reflected from the tooth surface. Prolonged exposure may cause retinal damage.
4. CAUTION: The Elipar FreeLight 2 unit generates high intensity light. The emitted light should be positioned directly above the material to be cured – prolonged light exposure of soft tissues (gingiva, oral mucosa, and skin) should be avoided as such exposure may cause damage or irritation. If applicable, cover such areas. If exposure of soft tissues cannot be avoided, adjust the polymerization process by shortening the polymerization times or by increasing the distance between the light guide exit and the material to be cured.
5. Condensation resulting from the unit being transferred from a cold to a warm environment may be a potential risk. Do not operate the unit until it is completely equilibrated to ambient temperature.
6. In order to avoid electric shock, do not introduce any objects into the unit with the exception of replacement parts handled in accordance with the Operating Instructions.
7. Use original 3M ESPE parts exclusively to replace defective parts in accordance with these Operating Instructions. The product guarantee does not cover any damage resulting from the use of third-party replacement parts.
8. Should you have any reason to suspect the safety of the unit to be compromised, the unit must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective unit. Safety may be compromised (e.g., if the unit malfunctions or is noticeably damaged).

Table of Contents	Page
Safety	1
Glossary of Symbols	2
Product Description	3
Fields of Application	3
Technical Data	3
Charger	3
Handpiece	3
Charger and Handpiece	3
Transport and Storage Conditions	3
Installation of the Unit	4
Factory Settings	4
Initial Steps	4
Charger	4
Light Guide/Handpiece	4
Insert the battery	4
Battery Charging	4
LED Display of the Charger	4
Table Holder for the Handpiece	4
Operation	5
Selection of Exposure Time	5
Activating and Deactivating the Light	5
Positioning the Light Guide	5
Removing and Inserting the Light Guide from/into the Handpiece	5
Measurement of Light Intensity	5
Low Battery Charge Display	6
Power-Down Mode	6
Acoustical Signals – Handpiece	6
Acoustical Signals – Charger	6
Operating Errors	6
Maintenance and Care	7
Replacement of the Battery	7
Handpiece/Battery Care	7
Cleaning the Light Guide	7
Cleaning the Charger, Handpiece, Table Holder and Glare Shield	7
Storage of the Handpiece during Extended Periods of Non-Use	8
Disposal	8
Customer Information	8
Limited Warranty – U.S. Only	8

9. Solvents, flammable liquids, and powerful heat sources must not be kept in the vicinity of the unit, as exposure to one of these factors may damage the plastic housing of the unit.
10. Do not allow cleaning agents to enter the unit, as this may short-circuit the unit or cause potentially dangerous malfunction.
11. Only factory-authorized technicians are authorized to open the housing of the unit and carry out any repairs.
12. Elipar FreeLight 2 must not be used in patients, with heart pacemaker implants or by users who have been advised to be cautious regarding their exposure to small electrical devices.
13. Do not use Elipar FreeLight 2 in patients with a history of photobiological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photosensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
14. Individuals with a history of cataract surgery may be particularly sensitive to the light exposure and should be discouraged from Elipar FreeLight 2 treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
15. Individuals with a history of retinal disease should seek advice from their ophthalmologist prior to operating the unit. In operating the Elipar FreeLight 2 unit, this group of individuals must take extreme care and comply with any and all safety precautions (including the use of suitable light-filtering safety goggles).
16. Prior to each use of the device, ensure that the emitted light intensity is sufficient to safely guarantee polymerization. Use the built-in light testing area for this purpose.

Glossary of Symbols


Attention, Consult Accompanying Documents

Type B Equipment

Class II Equipment – double insulated

Product Description

Elipar FreeLight 2, manufactured by 3M ESPE, is a high-performance light unit for intraoral polymerization of dental materials. The unit consists of a charger and a cordless handpiece powered by a rechargeable battery. The unit is designed for use on a table and cannot be wall-mounted.

The unit uses a high-performance Light Emitting Diode (LED) as the light source. In contrast to halogen light units, the unit emits light mainly in the wavelength range of 430 to 480 nm, e.g. the relevant range for camphor quinone-containing products.

The optimal match of the wavelength range to the intended purpose ensures that the polymerization performance is similar to that of halogen light-emitting units, albeit at lower light intensity. The polymerization performance of this light unit is so high, that the exposure times can be reduced by 50%, compared to a conventional halogen light unit.

Exposure time options: 5, 10, 15, or 20 seconds

The charger is equipped with an integral light intensity testing area.

The device is shipped with a turbo light guide with an 8 mm diameter light exit. It is not permissible to use the light guides of other units.

The "maxi fiber rod" with a diameter of 13 mm for larger areas (e.g. for fissure sealing), and the "proxi fiber rod" with a point-shaped light exit (e.g. for use in interproximal areas), are available as accessories.

Both the maxi fiber rod and the proxi fiber rod may only be used for the mentioned purposes and not for the polymerization of usual fillings, as otherwise complete polymerization can not be guaranteed.

The handpiece is equipped with a "power-down" function to minimize the unit's energy consumption. The handpiece switches to "power-down" mode once it is placed in the charger or if left unused for approximately 10 minutes outside the charger.

In stand-by mode, the charger consumes maximally 0.75 W. Starting in 2003, this value will be recommended for the stand-by mode by the EU according to the "Code of Conduct" on efficiency of external power supplies.

These Operation Instructions should be kept for the duration of product use.

Fields of Application

- Polymerization of light-curing dental materials with photo initiator for the wavelength range 430–480 nm

- Though the majority of light-curing dental materials are responsive in this range of wavelengths, you may wish to contact the manufacturer of the material in question to confirm the wavelength range.

Technical Data

Charger

Operating voltage: 120 V, 50/60 Hz

Power input: 0.1 A

Dimensions: depth: 210 mm
width: 95 mm
height: 60 mm

Weight: 555 g

Classification: Protection class II 

Handpiece

Power supply: Nickel-metal hydride battery, 4.8 V

Wavelength range: 430–480 nm

Light intensity: approx. 1000 mW/cm²

Duration of continuous use: min. 7 min (dependent on the ambient temperature before activation of temperature control)

Total operation time with new, fully charged battery: typically 20 minutes

Dimensions: diameter: 30 mm
length: 285 mm

Weight: 220 g

Charger and Handpiece

Time to charge empty battery: approx. 2 hours

Operating temperature: 10°C ... 40°C/59°F ... 104°F
Relative humidity: 30% ... 75%
Atmospheric pressure of 700 hPa to 1060 hPa

Total height with handpiece inserted in the charger: 190 mm

Transport and Storage Conditions:

- Ambient temperature range –20°C to +40°C/–4°F to +104°F

Subject to technical modification without prior notice.

Installation of the Unit

Factory Settings

The factory settings of the unit are 20 seconds operation time.

Initial Steps

Charger

- ▶ Please ensure first that the voltage stated on the rating plate corresponds to the existing main supply voltage. The rating plate is attached to the bottom of the unit.
- ▶ Place the charger on a level surface.
 - To protect the device from over-heating, do not obstruct the vent slots on the bottom of the unit.
- ▶ Connect the power cable of the charger to the power supply.
 - The green LED on the left side of the device is illuminated. This shows that the unit is ready for operation; please refer to the section, “LED display of the charger”.


Light Guide/Handpiece:

Never insert the handpiece in the charger unless the battery is inserted in the handpiece first!

- ▶ Steam Autoclave the light guide prior to first use.
- ▶ Then insert the light guide in the handpiece until it snaps into place.
- ▶ Place the enclosed glare shield on the light guide.

Insert the battery:

- ▶ Rotate the lid at the lower end of the handpiece counterclockwise until it hits the stop, and then remove the lid.


- ▶ To avoid malfunction, make sure the battery is properly inserted. Lay the handpiece upside down and insert the rechargeable battery in a fluent motion as directed by the arrow until you hear the battery snap into place.
- ▶ Replace the lid and rotate clockwise until it locks into place.

- ▶ In the event of malfunction, remove the rechargeable battery from the device and re-insert as described above.

Battery Charging:

- ▶ Prior to the first use, place the handpiece in the charger to fully load the new rechargeable battery.
 - The yellow LED of the charger illuminates after approximately 2 seconds; please see also section, “LED Display of the Charger”. The battery is fully charged once the yellow LED is turned off.
 - Note: Several cycles of charging and discharging may be necessary before the optimum cure time per charge is obtained.

LED Display of the Charger


LED green	LED yellow	Acoustical Signal		
On	Off	—	no handpiece inserted in the charger	→ charger ready for operation
On	Off	—	handpiece inserted in the charger	→ charging has been completed
On	On	—	handpiece inserted in the charger	→ rechargeable battery is being charged
On	Flashing	—	handpiece inserted in the charger	→ rechargeable battery is defective
On	Flashing	Ongoing beeps	handpiece inserted in the charger	→ charging contacts are wet

Table Holder for the Handpiece

While performing a procedure, the handpiece can be placed in a table holder.

Operation


Selection of Exposure Time

Exposure times of 5, 10, 15, and 20 seconds are available.

- ▶ For the compatible product at hand, use the exposure time as stated in the product's Instructions for Use and reduce to half.
- ▶ Select the exposure time by pressing the "sec" button.
 - The selected exposure time is indicated by the 4 green LEDs.
 - Each time the button is pressed, the setting advances to the next (higher) value; after 20 seconds the exposure time setting returns to 5 seconds. Keep the button depressed to scroll through the available settings.
 - The button for selection of the exposure time is deactivated during light activation.

Activating and Deactivating the Light

- ▶ Activate the light by briefly pressing the green Start button.
 - Initially, the LEDs show the preset exposure time: 4 illuminated LEDs represent 20 seconds of exposure time. After every 5 seconds of exposure, one of the LEDs is turned off (e.g., 3 LEDs correspond to 15 seconds of exposure time remaining, 2 LEDs to 10 seconds of exposure time remaining, etc.).
- ▶ If it is desired to turn off the light before the preset exposure time has completely elapsed, press the green Start button again.

Positioning the Light Guide

- ▶ Rotate the light guide into the optimal position for polymerization.
- ▶ To make full use of the light intensity provided, place the light guide as close to the filling as possible. Avoid directly contacting the filling material.
 - Keep the light guide clean at all times to obtain full light intensity.
 - **Damaged light guides substantially reduce the light intensity and must be replaced immediately. Sharp edges may cause serious injury.**

Removing and Inserting the Light Guide from/into the Handpiece

- ▶ To remove the light guide from the handpiece pull towards the front of the device.
- ▶ To put the light guide into the handpiece, push the light guide in until it snaps into position.


Measurement of Light Intensity

The light intensity can be reliably determined only with the charger of the Elipar FreeLight 2 unit. The light testing area is situated on the charger base. Testing the light intensity with any other unit produces erroneous results because of differences in the light sources used and arrangement of the components.

Caution: Measure the light intensity with the 8 mm turbo light guide only.

- ▶ If required, clean the test area with a wet cloth. Take care not to bend or damage the charging contacts.
- ▶ Without applying pressure, place the exit tip of the light guide level with the test area.
- ▶ Activate the lamp by pressing the green Start button.
 - All 5 blue LEDs come ON for approximately 1 second.
 - After this period of time, the number of illuminated LEDs is indicative of the measured light intensity: 5 LEDs = 100%, 4 LEDs = 80%, 3 LEDs = 60%, 2 LEDs = 40%, 1 LED = 20%.

- ▶ If the light intensity is below 80% (fewer than 4 LEDs are ON), check the light guide for contamination or defects.
- ▶ Proceed as follows:
 - 1. Clean the light guide if contaminated; please refer to the section on “Care”, or
 - 2. replace the light guide if defective, or
 - 3. if 1. and 2. fail to improve the situation, please contact 3M ESPE Customer Service or your retailer.

Low Battery Charge Display

After frequent use the battery charge of the device may drop to below approximately 10%, so that only a few more exposures can be performed without re-charging. Low battery charge is shown by several means:

- A short audible alarm is emitted 5x at the completion of an exposure, and repeated thereafter whenever any button is touched.
- The 5-second-LED of the handpiece begins to flash.
- ▶ The handpiece should be placed in the charger as soon as possible to re-charge the battery.

Power-Down Mode

Once the handpiece is placed in the charger, all internal functions and LEDs are automatically turned off as the handpiece switches to power-down mode. This reduces the power consumption of the rechargeable battery to a minimal level. Outside the charger, the handpiece is also switched to power-down mode if left unused for approximately 10 minutes.

- ▶ To terminate the power-down mode press one of the two buttons.
 - The power-down termination signal (two short audible signals) is emitted indicating that the

handpiece is ready for operation: the handpiece displays show the latest selected exposure mode and time settings.

Acoustical Signals – Handpiece

An acoustical signal is emitted

- every time a button is pressed,
- every time the light is turned ON or OFF,
- after 5 seconds of exposure time have elapsed (1 alarm signal), 2x after 10 seconds, 3x after 15 seconds.

Two acoustical signals are emitted

- every time the power-down mode is terminated by pressing any of the buttons.

A 2-second-error signal is emitted, if

- the temperature control is activated,
- the rechargeable battery is discharged to an extent that reliable polymerization is no longer ensured.


Acoustical Signals – Charger

Ongoing beeps are emitted when the handpiece is in the charger and the charge contacts are wet.

Operating Errors

Error	Cause ▶ Solution
The 5-second-LED flashes and a short signal is emitted 5x whenever a button is touched or the light is turned off.	The residual battery charge has dropped below 10%. ▶ Place the handpiece in the charger and re-charge the battery.
The ongoing exposure is interrupted (light off signal is emitted) followed by a 2 second-error signal; the device switches to power-down mode and resists further activation.	The battery lacks sufficient charge. ▶ Place the handpiece in the charger and re-charge the battery.
A 2 second-error signal is emitted upon pressing of the Start button.	The temperature control activated to protect the handpiece from overheating. The handpiece can be used again once it has cooled down. ▶ Allow the handpiece to cool down. <ul style="list-style-type: none"> - The light may be successfully activated once the handpiece has cooled down.

Error	Cause ▶ Solution
Ongoing beeps are emitted and the yellow LED is flashing fastly when the handpiece is put into the charger.	The charge contacts are wet. ▶ Dry the charge contacts. Do not bend the pins during drying.
The yellow LED of the charger flashes.	The battery is defective. ▶ Replace the battery.
The green LED of the charger fails to light up, even though the power cable is connected to a power outlet.	Power outlet carries no voltage. ▶ Use a different power outlet.
	Charger is defective. ▶ Have the charger repaired.
The light intensity is too low.	▶ Clean the light guide and the protecting glass in the light guide mounting hole (please refer to "Cleaning the Light Guide").

- ▶ Push the (replacement) battery into the handpiece proceeding along the direction shown by the arrow until the battery locks into place.
- ▶ Replace the lid and close by rotating clockwise.
- ▶ Insert the handpiece into the unit to completely charge the new battery for its first use.
 - The yellow LED of the charger illuminates after approximately 2 seconds; please see also "LED Display of the Charger". The battery charging process is complete once the yellow LED is turned off.
 - Note: Several cycles of charging and discharging may be necessary before the optimum cure time per charge is obtained.


Maintenance and Care

Replacement of the Battery

Never place the handpiece in the charger without the battery inserted in the handpiece.

Use 3M ESPE batteries only. The use of other manufacturer batteries or non-rechargeable or primary batteries is a potential hazard and may damage the device.

- ▶ Rotate the lid at the lower end of the handpiece counterclockwise until it hits the stop, and then remove the lid.
- ▶ Push a small spatula or similar aid between the cylinders of the battery until the snap-lock releases so that the battery protrudes slightly from the handpiece body.

- ▶ Remove the battery from the handpiece.

Handpiece/Battery Care

- ▶ Do not use any other chargers as this may damage the battery.
- ▶ The battery must not be immersed in water or thrown into open fires. Please see also "Safety, Section 5".

Cleaning the Light Guide

The light guide can be steam autoclaved. Do not sterilize by chemical means or hot-air.

- ▶ The light guide should be regularly wiped clean with a soft cloth. The spots of dried liquid should be wiped off the ends of the device especially before and after steam sterilization.
- ▶ Adhering polymerized composite should be removed with alcohol. A plastic spatula may help in removing the material.
 - To protect the surface of the device from scratching, do not use any sharp or pointed tools.
- ▶ The protecting glass in the light guide mounting hole can be cleaned with a dry cotton bud.

Cleaning the Charger, Handpiece, Table Holder and Glare Shield

- ▶ For disinfection of all components of the device, spray the disinfection agent onto a 4 x 4 gauze and proceed to disinfect the device. Birex SE, DisCide

Ultra and Sani-Cloth Plus are recommended surface disinfectants. Improper disinfection can cause material defects.

- **To avoid internal damage, do not spray disinfectant directly on handpiece or charger surface.** Disinfection agents must not enter the device. Apply surface disinfectant on the cool surface of the handpiece or charger.
- ▶ Dry residual disinfectants on the charger, the handpiece and the glare shield with a soft and fluff-free cloth, as they damage the plastic components.
- ▶ Clean the charger, the handpiece and the glare shield with a soft cloth and, if required, a mild cleaning agent.
 - Do not use solvents or abrasive cleaning agents as these may damage plastic parts of the device.
 - Cleaning agents must not enter the devices.
- ▶ Make sure that charge contact pins remain dry and are not contacted by metallic or greasy parts. Do not bend the charge contacts during drying. Wet charge contacts will cause an operating error (ongoing beeps and flashing yellow LED).

Storage of the Handpiece during Extended Periods of Non-Use

- ▶ If the handpiece is not to be used for an extended period of time (e. g., during vacation), fully charge the battery prior to departure or keep the handpiece inserted in the operational charger.
 - Otherwise, nearly discharged batteries may become excessively discharged due to the low degree of power consumption ongoing even in Power Down mode. This may damage the battery.
- ▶ Fully or nearly discharged batteries must be recharged as soon as possible.

Disposal

Your new device is equipped with a Nickel-metal hydride battery. This kind of battery is classified by the federal government as a non-hazardous waste. It is recyclable and safe for disposal in the normal municipal waste stream.

- ▶ Recycle or dispose of defective batteries and units in accordance with local legal regulations.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Limited Warranty – U.S. Only

3M ESPE guarantees your Elipar FreeLight 2 against defects in material and workmanship for two years from date of purchase. During the two-year warranty period, 3M ESPE will repair defects or replace a defective unit free if the unit was new when it was purchased, and if it is promptly returned to the following address:

3M Health Care Service Center
Suite 200, Bldg. 502-1W-01 Oakdale, MN 55128
1-800-292-6298

In Canada:
3M Canada Service Support Centre
80 Enterprise Drive South
London, Ontario N6N 1C2

Defects caused by misuse, neglect, accident, or abuse are not covered by this warranty.

This warranty does not cover lamps, fiber optic tips, and rechargeable batteries.

THIS WARRANTY IS EXPRESSLY IN LIEU OF AND EXCLUDES WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. 3M ESPE SHALL HAVE NO LIABILITY OF ANY KIND FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OR EXPENSES THAT ARISE DIRECTLY OR INDIRECTLY FROM THE USE OF THIS PRODUCT, OTHER THAN TO REPLACE OR REPAIR THE DEVICE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON. THE FOREGOING LIMITATIONS OF REMEDIES AND LIABILITIES MAY NOT BE CHANGED EXCEPT BY WRITTEN AGREEMENT SIGNED BY AN AUTHORIZED OFFICER OF 3M ESPE.

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3M, ESPE, Elipar and FreeLight are trademarks of 3M or 3M ESPE AG.
© 3M 2007. All rights reserved.
3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249.