3M ESPE	
Express™	
Vinyl Polysiloxane Impression Material	
Vinyl Polysiloxan Abformmasse	
Material de empreinte en vinyl polysiloxane	
Material e per impronte vinylpolisilossano	
Material de impresión de vinil polisiloxano	
Material para impressões de polisiloxano de vinilo	
Vinyl Polysiloxaan afdrukmaterial	
Αποτυπωτικό υλικό από βινυλοπολιλοξείδη	
Antytkrysmassa en vinylpolysiloxan	
Vinylly polysiloksiaani jäljennösmaeriaali	
Vinylpolysiloxaan atryksmateriale	
Vinylpolysiloksaan atryksmasse	
<div> <div>Instructions for Use</div> <div>Gebrauchsinformation</div> <div>Mode d'emploi</div> <div>Informazioni per l'uso</div> <div>Información de uso</div> <div>Instruções de uso</div> </div>	<div> <div>Грубітсінформаціе</div> <div>Οδηγίες χρήσεως</div> <div>Käytännösviite</div> <div>Κατάσταση</div> <div>Brugsanvisning</div> <div>Brugsanvisning</div> </div>
<div> <div>Made in Germany</div> <div>3M ESPE AG</div> <div>Dental Products</div> <div>8-2229 Seefeld - Germany</div> </div>	<div> <div>CE</div> </div>
<div> <div>3M ESPE</div> <div>Dental Products</div> <div>St. Paul, MN 55144-1000</div> </div>	<div> <div>3M ESPE</div> <div>Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249</div> </div>
<div> <div>3M ESPE, Express and Penta are trademarks of 3M or 3M ESPE AG.</div> <div>© 3M 2004. All rights reserved.</div> </div>	<div> <div>4400073/025/02</div> </div>

3M ESPE AG	
Dental Products	
8-2229 Seefeld - Germany	

3M ESPE	
Dental Products	
St. Paul, MN 55144-1000	
3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249	
3M ESPE, Express and Penta are trademarks of 3M or 3M ESPE AG.	
© 3M 2004. All rights reserved.	

ENGLISH	
Product Description	
Express® Penta™ H and Express® Penta™ H Quick are heavy-bodied consistency addition reaction silicone materials for use with the Pentamix™ mixing unit, all three products are manufactured by 3M ESPE. Express Penta H is hydrophobic and the Express Penta H and Express Penta H Quick foil bags are sealed with a PentaMatic sealing cap. The PentaMatic sealing cap automatically opens the foil bag once sufficient pressure is established by the plunger of the Pentamix.	
Express® Medium Body – Regular Set, Express® Light Body – Regular Set and Express® Light Body – Fast Set are medium-bodied and light-bodied consistency hydrophilic addition reaction silicone materials, manufactured by 3M ESPE, for direct intraoral application with the Garant Dispenser, manufactured for 3M ESPE.	

For details on all products mentioned below, please refer to the corresponding Instructions for Use. These Instructions for Use should not be discarded for the duration of product use.

Impression of inlay, onlay, crown, and bridge preparations.
Preparation Tray: Particularly suitable are rigid, non-perforated metal trays or custom plastic trays. Perforated trays may also be used. Blocking of the distal edge of the tray with a silicone or thermoplastic material provides for controlled flow properties affording beneficial effects especially in the distal molar area.
► For sufficient adhesion, apply a thin layer of VPS Tray Adhesive, manufactured by 3M ESPE, to the tray and allow to dry completely (at least 5 min., 15 min drying time is ideal).

Penta cartridge/poly bag:
► Place Express Penta H in the designated Express Penta H cartridge and Express Penta H Quick in the designated Express Penta H Quick cartridge only! Express Penta H and Express Penta H Quick poly bags must only be used according to the Fast Set, Medium Tip and Regular Set, respectively.
► Discard the first 3 cm of mixed paste extruded from newly filled cartridges prior to the first use for impression taking. The color of the dispensed paste must be uniform.

► If a new mixing tip was already installed when the cartridge was inserted, make sure the drive shaft is engaged prior to mixing.

Dispenser/cartridge:
► Place Express Medium Body – Regular Set, Express Light Body – Regular Set or Express Light Body – Fast Set in the Garant Dispenser.

Shortly before application, ensure that both cartridge openings are not clogged. When using a new cartridge, extrude a small amount of paste until both base and catalyst are fully mixed.

► Then attach a new Garant Mixing Tip – Yellow and, if necessary, attach an Intraoral Tip – Yellow.
► If desired, fill a separate Elastomer syringe directly from the mixing tip using an intraoral tip.

► Material that has set inside the mixing tip should not be extruded by force as this may damage the cartridge and mixing tip and result in the formation of leaks.

Retraction
Suitable retraction agents are e.g., aluminum hydroxide chloride and aluminum sulfate-based solutions.
► Keep the area dry from which an impression is to be taken.
► Remove the tooth or rings after using a subgingival retraction agent.
► Thoroughly remove the residue of the retraction agent by rinsing and drying.

Dosing and Mixing
► Express Penta H and Express Penta H Quick are automatically dosed and mixed in the Pentamix unit. The mixing ratio by volume is 5:1, base paste : catalyst.
► Express Medium Body and Light Body impression materials are dosed and mixed individually in the Garant Dispenser. The mixing ratio by volume is 1:1, base paste : catalyst.

Times

Processing time from start of mixing*	Intraoral setting time**
min.:sec	min.:sec
Express Penta H 1:30 (23°C/74°F)	5:00
Express Penta H Quick 1:00 (23°C/74°F)	3:30
Express Medium Body 1:30 (23°C/74°F)	5:00
Express Light Body 1:00 (23°C/74°F)	5:00
Express Light Body Regular Set 1:30 (37°C/95°F)	
Express Light Body Fast Set 1:00 (23°C/74°F)	3:30
	5:00 (37°C/95°F)

These are the processing times at 23°C/74°F or 37°C/95°F respectively. Higher temperatures shorten, and lower temperatures prolong, total processing times.

* Start of mixing = entry of paste into the mixing tip.
** For product combinations: the longer time stated is valid.

Impression Taking
Note
► At the onset of mixing, closely monitor the entry of paste into the mixing tip and ensure that the base paste and catalyst paste flow into the mixing tip (green and white paste with Express Penta H, light green and dark green paste with Express Penta H Quick, violet and light brown paste with Express Medium Body – Regular Set, green and blue paste with Express Light Body – Regular Set, and blue and grey paste with Express Light Body – Fast Set).
► The adhesive injection of the intraoral tip should be immersed in the prepared tooth may cause this material to set more rapidly than the tray material due to exposure to the higher temperature of the mouth. This may distort the impression.

One Step/Two Viscosities Technique
► Load the impression tray prepared with adhesive with Express Penta H or Express Penta H Quick.
► While the assistant loads the tray, the dentist may begin to apply Express Medium Body – Regular Set or Express Light Body – Regular Set over the preparation. Depending on the number of teeth, initiate application such that tray loading and application of material around the preparation are completed at the same time. During application, the tip of the intraoral tip should be immersed in the material at all times to avoid inclusion of air bubbles.

► Slowly place the loaded tray in the mouth parallel to the vertical axes of the prepared teeth, and take the impression without applying pressure. Avoid contact of the teeth with the tray.
► Once the setting time has elapsed, remove the impression from the mouth.

After Impression Taking
► Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene
► Place the impression in a standard disinfectant solution, e.g. Impreset™, manufactured for 3M ESPE, for the period of time recommended by the manufacturer, i.e. 10 min, if Impreset is used. Excessive disinfection may damage the impression.

** Impreset is not suitable for all countries.
► After disinfection, rinse the impression under running water for approx. 15 sec.

Pouring the Model
► To avoid introducing bubbles in the model, briefly pre-rinse with water and dry to use a silicone surfactant.
► Pour the impression into a commercial specialized stone plaster commencing after 30 min and no longer than 2 weeks. There are no other limitations on time.

► The impressions can be silver- or copper-coated.

Cleaning
► **Garant Dispenser:** Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135°C/275°F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
► **Impression Tray:** The adhesive can be removed from impression trays using acetone.

Notes
► The oxygen inhibition layer of composite materials, e.g. in fillings or post-core preparations, may impair the setting of silicone materials, and should therefore be removed completely before taking the impression.
► Keep the used mixing tip on top of the cartridge to serve as a seal until the next use. Removing and replacing the spent mixing tip may lead to carry-over of paste and the ensuing formation of hardened material.
► Express impression materials must not be combined with condensation reaction silicone materials as even traces of these materials impair setting.
► Disposable latex gloves impair the setting of impression materials. Hence, vinyl gloves are more suitable.

► The impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modeling.

Incompatibilities
► In suitable individuals, sensitization to the product can not be excluded. Use of the product should be discontinued and the product completely withdrawn from the patient, if sensitivity or allergic reactions are observed.

Storage and Stability
Store the product at 15–25°C/59–77°F. Do not use after the expiration date. Store impressions below 30°C/86°F.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your return of the product to 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Informational valid as of 06/04

One Step/Two Viscosities Technique
► Load the impression tray prepared with adhesive with Express Penta H or Express Penta H Quick.
► While the assistant loads the tray, the dentist may begin to apply Express Medium Body – Regular Set or Express Light Body – Regular Set over the preparation. Depending on the number of teeth, initiate application such that tray loading and application of material around the preparation are completed at the same time. During application, the tip of the intraoral tip should be immersed in the material at all times to avoid inclusion of air bubbles.

► Slowly place the loaded tray in the mouth parallel to the vertical axes of the prepared teeth, and take the impression without applying pressure. Avoid contact of the teeth with the tray.
► Once the setting time has elapsed, remove the impression from the mouth.

After Impression Taking
► Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene
► Place the impression in a standard disinfectant solution, e.g. Impreset™, manufactured for 3M ESPE, for the period of time recommended by the manufacturer, i.e. 10 min, if Impreset is used. Excessive disinfection may damage the impression.

** Impreset is not suitable for all countries.
► After disinfection, rinse the impression under running water for approx. 15 sec.

Pouring the Model
► To avoid introducing bubbles in the model, briefly pre-rinse with water and dry to use a silicone surfactant.
► Pour the impression into a commercial specialized stone plaster commencing after 30 min and no longer than 2 weeks. There are no other limitations on time.

► The impressions can be silver- or copper-coated.

Cleaning
► **Garant Dispenser:** Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135°C/275°F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
► **Impression Tray:** The adhesive can be removed from impression trays using acetone.

Notes
► The oxygen inhibition layer of composite materials, e.g. in fillings or post-core preparations, may impair the setting of silicone materials, and should therefore be removed completely before taking the impression.

► Keep the used mixing tip on top of the cartridge to serve as a seal until the next use. Removing and replacing the spent mixing tip may lead to carry-over of paste and the ensuing formation of hardened material.

► Express impression materials must not be combined with condensation reaction silicone materials as even traces of these materials impair setting.

► Disposable latex gloves impair the setting of impression materials. Hence, vinyl gloves are more suitable.

► The impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modeling.

Incompatibilities
► In suitable individuals, sensitization to the product can not be excluded. Use of the product should be discontinued and the product completely withdrawn from the patient, if sensitivity or allergic reactions are observed.

Storage and Stability
Store the product at 15–25°C/59–77°F. Do not use after the expiration date. Store impressions below 30°C/86°F.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your return of the product to 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Informational valid as of 06/04

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your return of the product to 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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