3M ESPE

Impregum™ Penta™

Polyether Impression Material

Polyether Abformmasse

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44000729/002

ENGLISH

Product Description
Impregum™ Penta™ is a medium-bodied consistency polyether impression material for the Pentamix™ mixing device, both products manufactured by 3M ESPE. The mixing ratio is 5 volumes base paste : 1 volume catalyst. The poly bags are sealed with a **PentaMatic™** sealing cap, manufactured by 3M ESPE. The **PentaMatic™** sealing cap automatically opens the poly bag once sufficient pressure is established by the plunger of the Pentamix. For details on Pentamix and accessories, Polyether Adhesive, Penta™ Elastomer syringe, Penta™ Mixing Tip, Permadyne™ Garant™ 2:1 and Impreset™, all products manufactured by or for 3M ESPE, please refer to the corresponding instructions for Use.

⚠ These instructions for use should not be discarded for the duration of product use.

Areas of Application

- Impression of inlay, onlay, crown, bridge, and veneer preparations
- Functional impressions
- Implant impressions

Preparation

- Block-out undercuts or areas where gingival recession exists to prevent the impression material from "locking" onto the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prosthesis. The block-out is especially important when using impression materials with high Shore hardness.

Impression trays

- Any impression trays generally used for precision impressions are suitable.
- For sufficient adhesion, apply a thin layer of Polyether Adhesive to the tray and allow to dry completely (at least 30–60 sec – 15 min drying time are optimal).

- Pentamix/Penta cartridge/poly bag:**
 - Place Impregum Penta in the designated cartridge only!
 - Poly bags equipped with a **PentMatic™** seal exclusively must be used with Pentamix™ Mixing Tips Red, manufactured by 3M ESPE.
 - Discard the first 3 cm of paste extruded from newly filled cartridges prior to the use for impression taking. The color of the dispensed paste must be uniform.
 - If a new mixing tip is installed when the cartridge is inserted, prior to mixing check whether the drive shaft is engaged.

Retraction
Solutions based on aluminum hydroxide chloride or aluminum sulphate are suitable retraction agents. Retraction with threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulphate or iron (III) sulphate may keep the setting of polyethers.
Areas from which impressions are to be taken should be kept dry.
► In subgingival preparations, hemostatic threads or rings may be used.
► Prior to taking the impression, completely remove all residual retraction agent by rinsing and drying.

Dosing and Mixing

- Dosing and mixing are performed automatically in the Pentamix.

Processing time from start of mixing*	Setting time from start of mixing*	Residence time in the mouth
min:sec	min:sec	min:sec
02:45	06:00	03:15

The working times shown are in conformity with ISO 4823:2000. Higher temperatures will shorten the working time, while lower temperatures will prolong it.

We recommend a working time of 01:45 min:sec and compliance with the indicated setting time from the start of mixing to insure an ideal impression result at all times under clinical conditions.

*Start of mixing = entry of paste into the mixing cannula.

Impression Taking

- If necessary for removal of the set impression from the mouth appropriately block-out deep undercuts in interdental areas.

Mono-phase technique

- For mixing around the preparation attach the Penta Elastomer syringe to the Penta Alloy Tip Rod of the Pentamix device, and fill the syringe.

► Then load the tray previously prepared with adhesive. The mixing tip should be kept immersed in the material at all times.

- Apply the material into the sulcus or cavity from below it, keeping the tip of the application nozzle immersed in the material and in contact with the tooth surface at all times.

- If a new mixing tip is used, the application of the material should be initiated only once the tray is half-loaded in order to insure that the tray can be positioned in the mouth immediately after insertion. Otherwise the syringe material may set more rapidly than the tray material causing distortion of the impression.

- Immediately after application around the preparation, position the loaded tray in the mouth, and hold it there without applying pressure until setting is complete.
- To remedy initial adhesion ("setting the valve"), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior direction. If this proves difficult, it may be necessary to carefully blow some air or water between the impression and the gingiva.

Functional Impression

- After preparing the custom tray with adhesive, load the tray with paste mix, position it in the mouth, and ask the patient to perform functional excursions.

- If required, use Permadyne Penta L or Permadyne Garant 2:1 for precise impression taking (correction).

After Impression Taking

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution commonly used for impressions, e.g. Impreset™, for the period of time recommended by the manufacturer, i.e., 10 min in the case of Impreset. Excessive disinfection may damage the impression.

- After disinfection, rinse the impression under running water for approx. 15 sec.

*Impreset is not available in all countries.

Model Preparation

- Prepare a cast from the impression with a specialized stone plaster no earlier than 30 min and no later than 14 days after impression taking.

- To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring. Do not use surfactants as these impair the quality of polyether impressions and do not use necessary!

- Polyether impressions can be silver-coated, whereas copper-coating is not possible.

Cleaning

- Paste that has not set may be removed with ethanol or by rinsing with water and soap.

- The adhesive can be removed from metal trays with acetone.

Notes

- At temperatures below 18°C/64°F, the viscosity of the pastes may increase sufficiently to make mixing in the device difficult. After keeping the pastes at 18°C/64°F for at least one day the processability is re-established without compromising quality.

- Direct exposure to sunlight and damp storage conditions may damage the impression.

- Polyether impressions should not be exposed to solvent-containing liquids, as this may result in swelling and impression distortion.
- 3M ESPE polyether materials may only be combined with 3M ESPE polyethers. The working and setting times of the individual products must be strictly observed.

Precautionary Measures

If an adverse reaction is observed in a patient, discontinue use of the product in this patient.

3M ESPE MSDSs can be obtained from www.mm.com or contact your local subsidiary.

Storage and Stability

Store the product at 18–25°C/64–77°F. **Do not refrigerate!** Store impressions dry for at least one day. Store impressions dry and below 30°C/86°F in the dark.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture for 30 days. NOT OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

3M ESPE does not warrant, by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, Garant, Impregum, Impreset, Penta, Pentamatic, Pentamix and Permadyne are trademarks of 3M or 3M ESPE AG.
© 3M 2007. All rights reserved.