## 3M ESPE

## Ketac™ Cem Aplicap™

### Glass Ionomer Luting Cement

### Glasionomer-Befestigungszement

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Made in Germany by	<b>CE</b>	<b>0123</b>
<b>3M ESPE AG</b>		
Dental Products		
D-82229 Seefeld - Germany		

**3M ESPE**
Dental Products
St. Paul, MN 55144-1000

3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada
1-800-634-2249.

3M, ESPE, Alkaliner, Aplicap, Capmix, Ketac and Rotomix are trademarks of 3M or 3M ESPE AG.

© 3M 2003. All rights reserved.

## ENGLISH

### Product Description

Ketac™ Cem Aplicap™, manufactured by 3M ESPE, is a shaded yellow fixing cement similar to dentin and comes in capsule form. The minimum amount capable of being dispensed from a capsule is 0.07 ml.

For details on Alkaliner™, Capmix™, Rotomix™ and Aplicap™ System, all products manufactured by or for 3M ESPE, please refer to corresponding instructions for Use. Instructions for Use should not be discarded for the duration of product use.

### Indication

Cementation of inlays, onlays, single crowns, posts and screws

Luting cement for use on orthodontic strips

Relinings

### Preparation

To provide for optimal adhesion, carefully clean enamel, dentin, and metal surfaces.

Blow the cavity dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry! **The cavity should be just so dry that the surface has a matt shiny appearance.** Excessive drying with compressed air can negatively affect adhesion or result in post-operative sensitivity after seating the restoration.

The metal surfaces should be dried well with compressed air. Avoid further contamination by saliva.

### Pulp Protection

The glass ionomer cement must not be applied directly to dentin – e.g. on inlays – situated close to the pulp or to the exposed pulp. Prior to taking the impression, cover areas close to the pulp using a hard-setting calcium hydroxide material, e.g. Alkaliner.

### Activation Mixing Application

Activate the capsules for 2 sec with the Aplicap Activator, mix using a high-speed mixer (approximately 4300 oscillations per min), e.g. Capmix or the Rotomix rotary mixer and apply with the Aplicap Applier. The Aplicap Activator and Applier are marked in red.

During application and the setting phase, water and saliva must be prevented from entering the working area. Apply a thin coat of cement to the inside of the crown and stump. Afterwards proceed with restoration. Avoid overfilling the crown.

Attention

Pauses between activation, mixing and application must definitely be avoided otherwise as soon as the material starts to set, removal from the capsule will be rendered difficult or even impossible. Hydrostatic pressure buildup in closely adjoining crowns may cause pulp discomfort.

### Times

The following time scale applies at a room temperature of 23°C/73°F:

	min <span> </span> : <span> </span> sec
Activation	0 <span> </span> : <span> </span> 02
Mixing in Rotomix	0 <span> </span> : <span> </span> 08
Mixing in high-speed mixer, e.g. Capmix	0 <span> </span> : <span> </span> 10
Working (from beginning of mix)	3 <span> </span> : <span> </span> 00
Setting (from beginning of mix)	7 <span> </span> : <span> </span> 00

The setting time is accelerated at room temperatures above 23°C/73°F. Exceeding the working time will cause loss of adhesion to enamel and dentin. Shorter mixing times will prolong the working and setting times, whereas longer mixing times will shorten them. On high-speed mixers the mixing time can be set from 8–15 sec and on the Rotomix from 7–12 sec.

### Removal of Excess

Remove excess using a Heidemann spatula and/or probe 7 min after beginning of mix.

### Incompatibilities

The product may cause allergic reactions in sensitive persons. If allergic reactions occur, use of the product should be discontinued and the product completely removed. Pulp complaints may be experienced in rare cases, in particular where working instructions are not properly observed.

### Storage and Shelf Life

Do not store the product above 25°C/77°F.

Do not use after the expiry date.

### Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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