<div> <div></div> <div>ESPE</div> </div> <p>Ketac™ Molar Aplicap™</p> <div> <div>Glass Ionomer Restorative Material</div> <div>Glasionomer Füllungsmaterial</div> <div>Matériau d'obturation en verre ionomère</div> <div>Materiale per otturazione vetro-ionomero</div> <div>Material de obturación de vidrio ionómero</div> <div>Material de obturacao de vidro ionómero</div> <div>Glasionomeer vulmateriaal</div> <div>Υαλοίονομερήs κoνία για εμφράξεις</div> <div>Glasjonomer fyllningsmaterial</div> <div>Lasi-ionomeeritäyteaine</div> <div>Glasionomer-fyldningsmateriale</div> <div>Glassionomer fyllingsmateriale</div> </div>

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Información per l'uso
Información de uso
Instruções de uso

<div> <div></div> <div>Made in Germany by</div> </div> <div> <div></div> <div>3M ESPE AG</div> </div> <div> <div></div> <div>Dental Products</div> </div> <div> <div></div> <div>D-82229 Seefeld - Germany</div> </div> <div> <div></div> <div>3M ESPE</div> </div> <div> <div></div> <div>Dental Products</div> </div> <div> <div></div> <div>St. Paul, MN 55144-1000</div> </div> <div> <div></div> <div>3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249.</div> </div> <div> <div></div> <div>3M, ESPE, Alkaliner, Aplicap, Capmix, Ketac and Rotomix are trademarks of 3M or 3M ESPE AG.</div> </div> <div> <div></div> <div>© 3M 2003. All rights reserved.</div> </div>	
<div> <div></div> <div>440007296561</div> </div> <div> <div></div> <div>CE</div> </div> <div> <div></div> <div>0123</div> </div>	

3M, ESPE, Alkaliner, Aplicap, Capmix, Ketac and Rotomix are trademarks of 3M or 3M ESPE AG.
© 3M 2003. All rights reserved.

<div> <div></div> <div>ENGLISH</div> </div>	
--	--

Product Description

Ketac™ Molar Aplicap™, manufactured by 3M ESPE, is a glass ionomer cement in a capsule. On account of its chemical adhesion to the dental enamel and dentin, a gentle preparation and an especially tight seal at the margins is possible. The product is available in the following shades: A1, A3, A4, and dark yellow opaque (DYO). The material can be applied without lining; it releases fluoride ions, and is radiopaque. The dispensable quantity of a capsule is at least 0.1 ml. For details on Ketac™ Conditioner, Alkaliner™, Capmix™, Rotomix™, Aplicap™ System and Ketac™ Glaze, all products manufactured by or for 3M ESPE™, please refer to the respective information for use. The information for use of a product has to be kept for the duration of the application.

Fields of application

Linings for uniplanar and multi-planar composite fillings
Core build-up prior to crown placement
Primary tooth fillings
Uniplanar fillings in non-occlusal load-bearing areas
Dental neck fillings, if aesthetics is not the prime consideration
Uniplanar and multi-planar temporary fillings

Preparation

Remove carious tooth substance only; undercuts are not required. Do not prepare thin bevelling margins. A wall thickness of 0.5 mm has to be secured for the material.

Pulp Protection

In order to prevent pulpal irritation cover areas in close proximity to the pulp with local application of a calcium hydroxide material (e.g. Alkaliner™, manufactured by 3M ESPE).

Conditioning

For optimal chemical bonding to the teeth the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner, manufactured by 3M ESPE, to the prepared surfaces and allowing the substance to react for 10 sec.

Then rinse with copious amounts of water.
Blow the cavity dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry!
The cavity should be just so dry that the surface has a matt shiny appearance. Excessive drying can result in post-operative sensitivity after filling.
Avoid further contamination.

Activation Mixing Application

Use the Aplicap activator to activate the capsules for 2 sec. Use a high-frequency mixer with abt. 4.300 oscillations/min, such as Capmix, or a rotation mixer Rotomix for mixing. Use the Aplicap applicler for application. Aplicap activator and applicler have been marked in red.

Attention
By all means prevent waiting times between activation, mixing, and application because the beginning of setting of the material will render difficult or prevent the application.

Ketac Molar Aplicap clings to metal instruments, and for this reason should be washed off with cold water before setting.

Times

The following times apply for an ambient temperature of 23°C/73°F:

	min : sec
Activation	0 : 02
Mixing in Rotomix	0 : 10
Mixing in high-frequency mixer e.g. Capmix	0 : 15
Working after start of mixing	2 : 00
Setting after start of mixing	4 : 30

Setting is accelerated over 23°C/73°F. An overextended working time will cause the loss of adhesion to the dental enamel and the dentin.
Shorter mixing times will extend, longer mixing times will shorten the working and setting times. The mixing times can be chosen between 8 and 15 sec in high-frequency mixers and between 7 and 12 sec in the Rotomix.

Filling protection

Glass ionomers go through two phases of setting, and are susceptible to water during the first setting period. Subsequently, dehydration may be possible, which can be avoided by applying Ketac Glaze directly after contouring. Apply Ketac Glaze to all free surfaces of the glass ionomer cement, and light-cure for 10 sec. The inhibition layer remaining after the setting of Ketac Glaze may obstruct the setting of the addition-type silicone impression materials. For this reason remove this layer or apply a polyether impression material.

Finishing

Further finishing with Arkansas stones, fine-grain diamonds, abrasive discs with successively decreasing grain size or silicone polishers should be started 4.5 minutes after start of mixing at the earliest. If desired, apply Ketac Glaze again for surface sealing and light-cure for 10 sec.

Incompatibility

The products mentioned may cause skin sensitisation in sensitive persons. If allergic reactions occur, use of the product should be discontinued and the product completely removed. In rare cases, especially in the event of non-adherence to the working instructions, pulplic complaints may occur.

Storage and shelf-life

Store the capsules in the blister packs.
Use within 1 month maximum once the blister pack has been opened!
Do not store the product above 25°C/77°F.
Do not use after the expiry date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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