3M ESPE	
Ramitec™ Penta™	
Polyether Bite Registration Material	
Polyether Bissregistriermaterial	
Ματériau à base de polyéther, pour empreintes occlusives	
Materiale per registrazioni occlusali a base di polietere	
Material para el registro oclusal a base de poliéter	
Material de poliéter para registros intramaxilares	
Polyether beetregistratiemateriaal	
Υλικό καταγραφής τήξεως από πολυαιθέρα	
Bettregistreringsmaterial av polyeter	
Polyeetteripohjainen -purennanrekisteröinti-materiaali	
Polyetherbidregistreringsmateriale	
Polyeter bittregistreringsmateriale	

Instructions for Use	Gebbruksinformatie
Gebrauchsinformation	Οδηγίες χρήσης
Mode d'emploi	Bruksanvisning
Informazioni per l'uso	Käyttöinfo
Informação de uso	Brugsanvisning
Instruções de uso	Bruksinformatjon

	
3M ESPE AG	CE
Dental Products	
D-82229 Seefeld - Germany	
3M ESPE	
Dental Products	
St. Paul, MN 55144-1000	
3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249	
3M, ESPE, Impresept, Penta, PentaMatic, Pentamix and Ramitec are trademarks of 3M or 3M ESPE AG.	
© 3M 2006. All rights reserved.	

	
---------------	---------------

ENGLISH

Product Description

Ramitec Penta is a polyether occlusal registration material for use with the Pentamix™ device, both products manufactured by 3M ESPE.The excellent flow properties of the product provide for exact occlusal reproduction without shifts. The material sets tough-elastically, is dimensionally stable, and resistant to heat and pressure. The foil bags are sealed with a **PentaMatic™** sealing cap, manufactured by 3M ESPE. The **PentaMatic** sealing cap automatically opens the bag once sufficient pressure is established by the plunger of the Pentamix. For details on all products mentioned below please refer to the corresponding Instructions for Use.

☞ Instructions for Use should not be discarded for the duration of product use.

Areas of Application

- Occlusal registration

Dosing and Mixing

► The material is dosed and mixed automatically in the Pentamix device.
► Foil bags equipped with a **PentaMatic** sealing cap must only be used with Penta Mixing Tips Red.
► Discard the first paste dispensed from a **new pair of plastic tubes** until the color of the paste emerging from the mixer is homogeneous.
► For filling of the Ramitec Application Syringe attach the rear opening of the syringe directly to the Penta Mixing Tip Red. The transparent body of the syringe permits the filling volume to be controlled. Filling is complete after approx. 20 sec.
► The used mixing tip may be kept on the cartridge to serve as a cap until the next use.

Application

► Apply a layer of Ramitec Penta with a thickness of approx. 5 mm to the row of teeth and have the patient close his/her mouth (final occlusion) or apply Ramitec Penta from the vestibular side in final occlusion.
For occlusal registration in the area of anterior teeth it is recommended that the paste is applied to the polished preparation and the opposing teeth.
► The impression can be removed 5 min after start of mixing. At this point the impression is of rubber-elastic consistency. The impression attains its final rigidity not until approx. 30 min after mixing.

Times	min:sec
Processing, from start of mixing	2:00
Setting, from start of mixing	5:00
Final rigidity, from start of mixing	30:00

Higiene

► Disinfect the occlusal impression using a standard disinfectant solution commonly used for the disinfection of impressions, e.g. Impresept™, for the period of time recommended by the manufacturer, i.e., 10 min in the case of Impresept. Excessive disinfection may damage the occlusal impression.
► After disinfection, rinse the occlusal impression under running water for approx. 15 sec.

* Impresept is not available in all countries.

Cleaning

► Clean the Ramitec Application Syringe only after the paste has completely set.
► Paste that has not set may be removed with ethanol or by rinsing with soap and water.

Notes

At low temperatures (below 18°C/67°F) the viscosity of Ramitec Penta may increase and make mixing in the device difficult. Hence, the foil bags should be stored at or above 18°C/67°F for one day prior to use. With this procedure a workable consistency of the paste is re-established without compromising quality.

Precautionary Measures

If an adverse reaction is observed in a patient, discontinue use of the product in this patient.

3M ESPE MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

Storage and Stability

Store the product at 18–25°C/64–77°F. **Do not refrigerate!** Do not use after the expiry date. Store occlusal impressions below 30°C/86°F in the dark.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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