ENGLISH

DESCRIPTION

Scotchbond™ Etchant delivery system, manufactured by 3M ESPE, is comprised of a syringe having a luer lock end and interlocking, disposable dispensing tips. The syringe contains 35 weight percent phosphoric acid etchant. The system is indicated for intraoral delivery of the acid etchant for enamel etching in the acid etch technique.

PRECAUTIONS FOR DENTAL PERSONNEL AND PATIENTS

Etchant Precautions: Etchant contains 35 weight % phosphoric acid. Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes and skin. If accidental contact occurs, flush immediately with large amounts of water. For eye contact, also consult a physician.

See Material Safety Data Sheet (MSDS) for additional cautionary information. You may obtain the current MSDS by visiting the website: <http://www.3m.com/MSDS> or contacting your 3M ESPE Dental Products representative.

TECHNIQUE

1. Protective eyewear for patients and dental staff is recommended when using the delivery system.
2. **Prepare Delivery System:** Remove clear cap from etchant syringe (Figure 1) and **SAVE**. Twist a blue disposable dispensing tip securely onto the syringe (Figure 2). Holding the tip away from the patient and any dental staff, express a small amount of etchant onto a dispensing pad or a 2 x 2 gauze to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and again express etchant. If clog remains, discard dispensing tip and replace with a new one. Bend the dispensing tip to a desired angle. Place bend midway along tip (Figure 3). Do not bend dispensing tip at its hub.
3. **Clean Enamel:** If enamel to be etched is not freshly cut, clean with a pumice-water slurry. Rinse well. Prophy pastes containing oil or fluoride should not be used as these additives may interfere with etching.
4. **Isolate:** A dry operating field is essential for the acid etch procedure. Rubber dam is the preferred method of isolation.
5. **Enamel Etching:** Using the delivery system, gently extrude etchant directly onto enamel surfaces to be etched (Figure 4). Leave etchant in place for 15 seconds.
6. **Rinse Etched Surfaces:** Thoroughly rinse for 15 seconds. Remove rinse water with suction. Do not allow patient to rinse. If saliva contacts the etched surfaces, re-etch for five seconds and rinse. Re-isolate the etched surfaces if using cotton rolls.
7. **Dry Etched Enamel:** Thoroughly dry the etched surfaces. Air should be oil and water free. The dry etched surfaces should appear frosty white. If not, repeat steps three, four, five and six.

DO NOT ALLOW THE ETCHED SURFACES TO BE CONTAMINATED BY ANYTHING. If contamination occurs, re-etch, rinse and dry as above.

8. Immediately proceed with the bonding procedure.

STORAGE & USE

1. The etchant is designed for use at room temperatures of approximately 70 - 75° F (21-24° C).
2. Shelf life at room temperature is 36 months. See outer package for expiry date.
3. **Delivery System Storage:** Remove used dispensing tip (Figure 5) and discard. Twist on clear storage cap (Figure 6). **Storage of the delivery syringe with a used dispensing tip or without the storage cap will allow drying of the etchant and consequent clogging of the system.** Replace storage cap with a new dispensing tip at next use.
4. If desired, the etchant may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.
5. If a liquid etchant is desired, the Etchant may be dispensed into a dappen dish and stirred to increase its fluidity.
6. **Disinfection:** Discard used dispensing tip. Replace syringe storage cap. Disinfect the capped syringe in the same manner as nonimmersible handpieces, air/water syringe and ultrasonic scalers following American Dental Association (ADA) and Centers for Disease Control (CDC) recommendations.¹

¹Council on Dental Materials, Instruments, and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, and Scotchbond are trademarks of 3M or 3M ESPE AG. © 3M 2002. All rights reserved. 44-0023-3962-8-A
3M ESPE AG