<b>3M ESPE</b>
<b>Visio™ Bond</b>
<b>Light-curing Bonding-Agent</b>
<b>Lichthärtendes Bonding-Material</b>
<b>Agent d’adhérence photopolymérisable</b>
<b>Materiali de adesione fotoindurente</b>
<b>Agente aderente fotopolimerizabile</b>
<b>«Bonding» Fotopolimerizável</b>
<b>Lichtuithardend adhesief</b>
<b>Φωτοπολυμεριζόμενος συγκολλητικός παράγοντας</b>
<b>Ljushärdande bondingvätska</b>
<b>Valokovetteinen sidosaine</b>
<b>Lyshærdende adhævis</b>
<b>Lysherdende kontaktformidler</b>

Instructions for Use	Gebbruksinformatie
Gebrauchsinformation	Πληροφορίες χρήσης
Mode d’emploi	Bruksanvisning
Informazioni per l’uso	Käyttöohjeet
Información de uso	Bruksanvisning
Instruções de uso	Bruksveiledning

**3M ESPE AG**
Dental Products
D-82229 Seefeld – Germany

**3M ESPE**
Dental Products
St. Paul, MN 55144-1000

3M ESPE Technical Hotline/MSDS Information in U.S.A. and Canada 1-800-634-2249

3M, ESPE, CoJet, Elipar, Sinfony and Visio are trademarks of 3M or 3M ESPE AG.
© 3M 2002. All rights reserved.

ENGLISH
---------

**Product Description**

Visio Bond, manufactured by 3M ESPE, is a one-component bonding material for intraoral repair of ceramic and composite restorations after prerequisite pre-treatment with CoJet™ Prep Microblaster and CoJet™ Sand Blast-Coating Agent, and ESPE™ Sil Silane Solution, all products manufactured by or for 3M ESPE. For details on all products mentioned below please refer to the corresponding Instructions for Use.

☞ Instructions for Use should not be discarded for the duration of product use.

**Areas of Application**

- Bonding of ceramic or composite restorations to a composite filling material in intraoral repair applications.

**Application**

- Visio Bond also polymerizes if exposed to daylight or artificial light and therefore must be protected from exposure to strong light sources.
- Use the CoJet Prep Microblaster and CoJet Sand Blast-coating agent to pretreat the site of the restoration needing to be repaired.
- Then apply ESPE Sil silane solution and allow to dry for 30 sec intraorally.
- Apply Sinfony™ Opaquer, manufactured by 3M ESPE, to exposed metal and light-cure with a halogen light unit.
- Dose Visio Bond onto a mixing pad directly prior to processing. Immediately reseal the bottle.
- Use a disposable brush to apply a thin layer of the bonding agent to the silanized ceramic or silanized composite material.
  - Freshly applied and light-cured opaquer does not need to be wetted, though wetting is not detrimental.
- Light-cure the bonding material for 20 sec.
- Subsequently, apply the composite filling material in accordance with the manufacturer’s recommendations.

**Incompatibilities**

In susceptible individuals, sensitization to the product cannot be excluded. If allergic reactions occur, use of the product should be discontinued and the product completely removed.

**Storage and Shelf-life**

Do not store the product above 25°C/77°F. Do not use after the expiration date.

**Customer Information**

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE’s sole obligation shall be repair or replacement of the 3M ESPE product.

**Limitation of Liability**

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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