## General

Z100™ Restorative material, manufactured by 3M ESPE, is a visible-light activated, radiopaque, restorative composite. It is designed for use in both anterior and posterior restorations. The filler in Z100 Restorative is zirconia/silica. The inorganic filler loading is 66% by volume with a particle size range of 3.5 to 0.01 micron. Z100 Restorative contains BIS-GMA and TEGDMA resins. Adper™ Scotchbond™ adhesive, manufactured by 3M ESPE, is used to permanently bond the restoration to the tooth structure. The restorative is available in a variety of shades. It is packaged in traditional syringes and single dose capsules.

**Precautions For Patients and Dental Personnel**

**1. ETCHANT PRECAUTIONS:** Contact with eyes or skin may cause burns. In case of contact, flush area with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection.

**2. COMPOSITE PASTE PRECAUTIONS:****Precautionary Information for Patients:**

Avoid use of this product in patients with known acrylate allergies. This product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

**Precautionary Information for Dental Personnel:**

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M ESPE MSDSs can be obtained from [www.mmm.com](http://www.mmm.com) or contact your local 3M subsidiary.

**INSTRUCTIONS:** The following instructions are separated into three sections: preliminary, anterior placements and posterior placements.

**I. PRELIMINARY**

**A. PROPHY:** Teeth should be cleaned with pumice and water to remove surface stains.

**B. SHADE SELECTION:** Before isolating the tooth, select appropriate shade of material. Some hints for choosing the correct shade are listed below.


**SHADE:** Teeth are not monochromatic. Each of the three areas of the tooth has a characteristic color.

**GINGIVAL AREA:** If the restoration is in the gingival area of the tooth, make note of the amount of yellow color it contains.

**BODY AREA:** Observe the body of the tooth and note the characteristic color it contains- is it gray, yellow, or universal?

**INCISAL AREA:** Study the incisal edges of the tooth and the neighboring teeth. Do the incisal edges have a blue or gray color to them? Note how far the translucency extends. This should be duplicated in the repaired tooth.

**THINNER RESTORATIONS:** The amount of color that a restoration will have is partly due to the thickness of the restoration. Thus, if a shade match is taken from a thick

**NOTE:** The matrix may be placed following the enamel etching and adhesive application steps if preferred.

**D. ETCHING:** Etch according to the instructions of the etchant that is used.

**E. PRIMING:** Prime according to the instructions of the Adper Scotchbond adhesive product that is used.

**F. ADHESIVE APPLICATION:** Apply the adhesive according to the directions of the Adper Scotchbond adhesive product that is used.

**G. ADHESIVE CURING:** Cure the adhesive coating the appropriate time according to the Adper Scotchbond adhesive product used. Expose its entire area to visible light from a 3M ESPE light, manufactured by 3M ESPE, or other dental visible light curing unit of comparable intensity.

**H. DISPENSING THE COMPOSITE:** Follow the directions corresponding to the dispensing system chosen.

**1. SYRINGE**

a. Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of the restorative when dispensing is completed, turn the handle counterclockwise a half turn to stop paste flow. Immediately replace syringe cap. If not used immediately, the dispensed material should be protected from light.

b. Place restorative into the cavity using a nonmetallic placement instrument. Overfill the cavity to permit extension of composite beyond cavity margins. Contour and shape with appropriate composite instruments.

**2. SINGLE DOSE CAPSULE**

(See separate 3M ESPE Delivery System, manufactured for 3M ESPE, Instructions)

**I. CURING:** Cure with a 3M ESPE light or other dental visible light curing unit of comparable intensity. Hold the light exit tip as close to the restorative material as possible. A cure time for each shade is listed below.

| SHADE | THICKNESS (mm) | TIME (sec) |
|-------|----------------|------------|
| A1.0  | 2.5            | 40         |
| A2.0  | 2.5            | 40         |
| A3.0  | 2.5            | 40         |
| A3.5  | 2.5            | 40         |
| B2.0  | 2.5            | 40         |
| B3.0  | 2.5            | 40         |
| C2.0  | 2.5            | 40         |
| C4.0  | 2.0            | 40         |
| D3.0  | 2.5            | 40         |
| P     | 2.5            | 40         |
| I     | 2.5            | 40         |
| A4.0  | 2.0            | 40         |
| CY    | 2.0            | 40         |
| CG    | 2.0            | 40         |
| UD    | 2.0            | 40         |

**J. FINISHING:** Contour restoration surfaces with fine finishing diamonds, burs, or stones. Contour proximal surfaces with Sof-Lex™ Finishing Strips, manufactured for 3M ESPE.

**K. ADJUST OCCLUSION:** Check occlusion with a thin articulating paper. Centric and lateral excursion contacts should be examined. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

**L. POLISHING:** Polish with Sof-Lex Discs, manufactured by 3M ESPE, and strips and with white stones or rubber points where discs are not suitable.

**III. POSTERIOR RESTORATIONS**

**A. CAVITY PREPARATION:** Cavity design requirements are essentially a conventional preparation with refinement of the cavosurface margin for enhancement of acid etching. No residual amalgam or other base material should be left in the internal forms of the preparation which would interfere with light transmission and the hardening of the restorative.

**B. PULP PROTECTION:** Use Vitrebond Liner/Base to base areas of deep cavity excavation. If a pulp exposure has occurred, use a minimum amount of calcium hydroxide followed by an application of Vitrebond Liner/Base. In deep restorations without pulp exposure, only the Vitrebond Liner/Base is needed. See Vitrebond instructions for details.

**C. PLACEMENT OF MATRIX:** Place a thin dead-soft matrix band and insert wedges firmly. Burnish the matrix band to establish proximal contour and contact area. Adapt the band to seal the gingival area and avoid overhang.

**D. ETCHING:** Etch according to the instructions of the etchant that is used.

**E. PRIMING:** Prime according to the instructions of the Adper Scotchbond adhesive product that is used.

**F. ADHESIVE APPLICATION:** Apply the adhesive according to the instructions of the Adper Scotchbond adhesive product that is used.

**G. ADHESIVE CURING:** Cure according to the instructions of the Adper Scotchbond adhesive that is used.

**H. DISPENSING THE COMPOSITE:** Follow the directions corresponding to the dispensing system chosen.

a. Incremental placement of the material is recommended. Avoid intense light in the working field. Using a plastic placement instrument or the capsule system, place the first increment of material not to exceed 1.5mm in thickness into the proximal box preparation. Using a condensing instrument, adapt the material to the internal surfaces of the box preparation and to the matrix and cure. Complete the box preparation with incremental placements not exceeding 2.5mm in thickness. Cure each increment separately.

b. Fill the remainder of the cavity preparation in layers no greater than 2.5mm. For large cavity preparation, it is recommended that the final occlusal placement be done in vertical increments buccal to lingual. Contour to provide proper contact area and occlusal anatomy.

**I. FINISHING:** Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with Sof-Lex Finishing Strips.

**J. ADJUST OCCLUSION:** Check occlusion with thin articulating paper. Centric and lateral excursion contacts should be examined. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

**K. POLISHING:** Polish with Sof-Lex discs and strips and with white stones or rubber points where discs are not suitable.

**IV. Z100 RESTORATIVE INDIRECT PROCEDURE FOR INLAYS****Dental Operatory Procedure**

**1. Shade selection:** Choose the appropriate shade(s) of Z100 Restorative prior to isolation. If the restoration is of sufficient depth, use of the opaque shade is recommended. Use of Incisal shade on the occlusal surface will help to achieve aesthetic appearance.

**2. Preparation:** No changes are required in the traditional inlay preparation technique. However, the preparation should allow for a minimum of 2mm depth and 2mm width of Z100 Restorative.

**3. Impressioning:** After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. Either the Express™ Impressioning System or the Imprint™ Impression Material, manufactured by 3M ESPE, may be used.

**LABORATORY PROCEDURE**

1. Pour the impression of the prepared inlay with die stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used.

2. Separate the cast from the impression after 45 to 60 minutes. Pindex, or Di-Lok,

44-0007-4343-3-B

Dental Products

3M ESPE

10/05

Instructions for use  
Mode d'emploi  
Istruzioni per l'uso**3M ESPE****Z100™ Restorative**

Z100™ Restorative

Z100™ Matériau de restauration

Materiale da Restauo Z100™

shade guide tooth for a thin restoration, an incorrect shade may be chosen. For a thin restoration, the wedge end of the shade guide should be used.

**MOCK UP:** Place the material, in the chosen shade on the unetched tooth. Manipulate the material to approximate the thickness and site of the restoration. Cure. Have several people evaluate the mock-up under different light sources. Flip the restorative off the unetched tooth with an explorer. If the shade did not match, choose another and repeat the mock-up. If the shade was acceptable, continue the tooth preparation with isolation.

**C. ISOLATION:** A rubber dam is the preferred method of isolation. Cotton rolls plus an evacuator can also be used.

**II. ANTERIOR RESTORATIONS**

**A. CAVITY PREPARATION:** Use conventional cavity preparations for all Class III, IV and Class V restorations.

**B. PULP PROTECTION:** Use Vitrebond™ Light Cure Glass Ionomer Liner/Base, manufactured by 3M ESPE, to base areas of deep cavity excavation. If a pulp exposure has occurred, use a minimum amount of calcium hydroxide followed by an application of Vitrebond Liner/Base. In deep restorations without pulp exposure, only Vitrebond Liner/Base is needed. See Vitrebond instructions for details.

**C. PLACEMENT OF MATRIX:** Mylar strips and crown forms may be used to minimize the amount of excess material used.

and base the cast as you would for a typical crown and bridge procedure. Mount or articulate the cast to its counter model to an adequate articulator.

3. If a second impression was not sent, pour the inlay impression using the same impression registration to be used as a working cast.

4. Section out the preparation with a laboratory saw and trim away excess, or expose so that the margins can be easily worked. Mark the margins with red pencil if needed. Add a spacer at this time if one is being used.

5. Soak the die in water, then with a brush, apply a very thin coat of separating medium (Liquid Foil) to the preparation, let it dry somewhat, then add another thin layer.

6. Add the first third of the resin to the floor of the preparation, stay short of the margins, light cure for 40 seconds.

7. Add second third of resin, capturing the medial/distal lower two thirds of the margin. Allow for the last third (incisal) to include the contact areas, light cure for 40 seconds.

8. Place the die back into the articulated arch add the last third of incisal resin to the occlusal. Over fill slightly mesially/distally, and occlusally. This will allow for the m/d contacts and for the proper occlusal contact when the opposing is brought down into occlusion with the still uncured incisal resin. Light cure for only ten seconds, then take the die out so it does not stick to the adjacent dies. Finish the curing process.

9. With the occlusal contacts already established, begin removing the excess resin material from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.

10. Removing the inlay from the stone die must be done with care. Break off small amounts of the die from around the restoration. The die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

11. Using master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, then polish.

#### Dental Operatory Procedure

1. Roughen the interior surfaces of the inlay.

2. Clean the Z100 Restorative inlay in a soap solution in an ultrasonic bath.

3. Cementation: Cement the inlay using a resin cement, e.g., RelyX™ ARC adhesive resin cement, and Adper Scotchbond Multi-Purpose Plus dental adhesive system, manufactured by 3M ESPE, by following manufacturer's instructions.

#### ADDITIONAL NOTES

1. Sensitivity – some patients may experience transitory postoperative sensitivity. The risk of sensitivity can be minimized by the following measures:

a. Remove minimal tooth structure.

b. Use proper isolation. Use of a rubber dam is highly recommended.

c. Adequate pulp protection. Use a glass ionomer cavity liner/base on appropriate dentinal surfaces.

d. Place restorative material in increments, curing each increment separately.

e. Adequately cure restorative according to instructions for shade and thickness of restorative and light exposure time.

f. Adjust occlusion carefully. Check for hyperocclusion, particularly in lateral excursion contacts.

#### STORAGE AND USE

1. Do not expose restorative materials to elevated temperatures or intense light.

2. Do not store materials in proximity to eugenol-containing products.

3. The composite pastes are designed for use at room temperature of approximately 21-24°C or 70-75°F. Shelf life at room temperature is 36 months for syringes; 36 months for capsules. See outer package for expiry date.

4. DO NOT remove capsules from the bottle until ready to use paste.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

#### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, Adper, Scotchbond, Z100, Vitrebond, RelyX, Sof-Lex, Express and Imprint are trademarks of 3M or 3M ESPE AG.
© 3M 2005. All rights reserved.