

# Filtek™ Z250 Universal Restorative

## Ordering Information

ordering

Shades	Item Numbers	
	Syringe	Capsule
A1	6020A1	6021A1
A2	6020A2	6021A2
A3	6020A3	6021A3
A3.5	6020A3.5	6021A3.5
A4	6020A4	6021A4
B0.5	6020B0.5	6021B0.5
B1	6020B1	6021B1
B2	6020B2	6021B2
B3	6020B3	6021B3
C2	6020C2	6021C2
C3	6020C3	6021C3
C4	6020C4	6021C4
D3	6020D3	6021D3
Incisal	6020I	6021I
Universal Dentin	6020UD	6021UD

3M, ESPE, Filtek and Silux Plus are trademarks of 3M or 3M ESPE AG. Charisma is a trademark of Heraeus Kulzer, Inc. Caulk, Dispersalloy, Esthet-X, TPH, TPH Spectrum, and SureFil are trademarks of Dentsply. Herculite XRV and Prodigy are trademarks of Kerr Corporation. Vitalescence is a trademark of Ultradent Products, Inc. Tetric Ceram is a trademark of Ivoclar Vivadent.

3M ESPE Technical Hotline: 1-800-634-2249

3M ESPE Web site: [www.3MESPE.com](http://www.3MESPE.com)

**3M ESPE**

### Dental Products

3M Center  
Building 275-2SE-03  
St. Paul, MN 55144-1000  
USA

### Dental Products

#### 3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1 800 265-1840 ext. 6229



Minimum 10%  
Post-Consumer Fiber

Printed in U.S.A.  
© 3M 1999-2002  
70-2009-2189-1  
44-0023-3777-0

# Filtek™ Z250 Universal Restorative

Reduced shrinkage, proven performance  
and excellent esthetics.



### KEY FEATURES

- LOW SHRINKAGE
- EXCELLENT ESTHETICS
- ALL-AROUND VERSATILITY



**3M ESPE**

# Filtek™ Z250 Universal Restorative



6020 - Filtek™ Z250  
Universal Restorative  
System Syringe  
Introductory Kit



6021 - Filtek™ Z250  
Universal Restorative  
System Capsule  
Introductory Kit

The lowest polymerization shrinkage, as shown in Chart 1, is only the beginning. With its high performance, strength, and easy handling, this micro-hybrid restorative fills a variety of needs—with beautiful results.

Filtek Z250 Universal Restorative offers:

**Lowest Polymerization Shrinkage** (See Chart 1)—Less stress on the tooth, less sensitivity potential for the patient.

**Speed and Convenience**—Easy-to-use, direct-delivery capsule. Cures in just 20 seconds\* and finishes to a smooth surface.

**Easy Handling**—Exceptional resistance to slumping.


**Outstanding Versatility**—Useful for any area of the mouth.

**Excellent Esthetics**—Available in 15 shades to meet your needs.

\*Shades UD, C4, and B0.5 require a 30-second cure for a 2mm increment.


## 1. Volumetric Shrinkage

*Lowest shrinkage of major competitors. 3M ESPE R&D conducted a shrinkage test on a variety of market-leading restoratives. Of those tested, Filtek™ Z250 Universal Restorative demonstrated the least amount of polymerization shrinkage.*



## 2. Wear Resistance

*Excellent strength and durability. Filtek™ Z250 Universal Restorative demonstrated high wear resistance when compared to other market-leading restoratives.*



Source: Dentsply International Esthet-X Brochure, 1999.  
[www.caulk.com/Profiles/EsthetXBrochure.htm](http://www.caulk.com/Profiles/EsthetXBrochure.htm)