

## Vitremer™ Core Buildup/Restorative System

### Dispensing:

- Dispense an equal number of powder scoops and liquid drops.

### Dispensing the powder:

- Fluff the powder.
- Dispense level scoops.
- Dispense equal numbers of powder scoops and liquid drops.



1

### Dispensing the liquid:

- Keep the tip of the vial clean.
- Hold the vial in a vertical orientation and squeeze.
- Avoid contacting the mixing pad with the tip.



2

### Mixing:

- Use a large cement spatula.
- First mix one-half of the powder into the liquid, then add the remaining powder.
- Mix aggressively.
- Continue mixing until all of the powder is incorporated into the liquid.
- Time the mixing. Do not mix longer than 45 seconds.
- Vitremer is designed to have a thick consistency.



3