

PLANNING A FIXED HYBRID RESTORATION

Pre-Surgical Appointment (Initial Steps identical for any Implant Restoration)

- ◆ Pano/CBCT Scan
- ◆ Case Planning Meeting
 - Surgical Dr., Restorative Dr., Laboratory Team Meeting
 - Determine current Restorative Interocclusal Space available
 - Review Proposed Solutions and Determine Projected Outcome/ Concerns
- ◆ Set-up or Wax-up Proposed Solution
 - Initial Impressions, Bite Registration, Mounting
- ◆ Approved Final Restoration Wax-up



Zirconia

- ◆ Scan Wax-up and CAD Design Restoration
 - Mill (CAM), Infiltrate, Sinter, and Finish Restoration
 - If Porcelain bonded to Zirconia Restoration, "Green State" cutback prior to Infiltration, Sinter and Build Up Porcelain

Denture

- ◆ Conversion
 - Existing Denture (upper/lower/both upper and lower) or impressions for an immediate Denture(s)
 - If required Immediate Dentures will be constructed
- ◆ Bite Registration- Kois Facial Analyzer/Gothic Arch Tracing/Bite Rims contoured with Midline, High Lip line, Corners (Cuspids)
- ◆ Shade
- ◆ Photo's/Video (Full Face, Relaxed Lips, Phonetic Sounds (F's and S's), Big Smile, Full Face with retractors, Side Profile) Can be done at the lab.

RVDA Removable Hybrid Pricing

Standard Hybrid AO4	\$3500
Hybrid w/Conversion & Interim Denture	\$4500
Fixed Removable	\$6500
Zirconia	\$5500

RVDA
ROGUE VALLEY DENTAL ARTS[®]

(541) 779 - 6503